

Compatibility Chart Expell Sealing Membranes with PCR/qPCR 96 Well Plate

CAPP Cat. No.	5100219C	5100220C	5100221C	5100319C	5100318C	5100317C
Image						
Standard Well (0.2 mL, height 21 mm)	■	■	■			
Low Profile Well (0.1 mL, height 16 mm)				■	■	■

CAPP Sealing Membranes			PCR/qPCR 96 Well Plate						
Cat. No.	Description	Image	CAPP Cat. No.	5100219C	5100220C	5100221C	5100319C	5100318C	5100317C
5100400C	Expell PCR 96-well plate Optical Sealing Membrane, ultra clear, pressure-sensitive adhesive film, operating temperature -40-120°C, Real-Time PCR (qPCR) use			■	■	■	■	■	■
5100401C	Expell PCR 96-well plate Sealing Membrane, strong adhesive, operating temperature -40-120°C, PCR use			■	■	■	■	■	■
5100402C	Expell qPCR 96-well plate Optical Sealing Membrane, ultra clear, strong adhesive, operating temperature -40-120°C, Real-Time PCR (qPCR) use			■	■	■	■	■	■
5100403C	Expell PCR 96-well plate Aluminium Membrane, strong adhesive, operating temperature -80-120°C, PCR and storage use, piercable			■	■	■	■	■	■