

CAPP *Rhythm*

MECHANICAL REPEATER

Features:

- Reliable ergonomics for your lab routine without sacrificing ultimate precision
 - Covers the complete set of dispensing volumes 1-5000 μ L
 - Single-handed volume selection, loading and dispensing
 - Minimum maintenance requirements
 - Robust design and construction
 - Almost no wear on parts
 - Lightweight – only 105g
 - One year warranty



Ordering Information

Cat. No.	Description
PR-50	CAPPRhythm Mechanical Repeater, 0.1-50mL

CAPP *Rhythm* Dispensing Volumes

Setting	1	2	3	4	5
TIPS/ Steps	48	23	15	11	8
0.05mL	1 μ L	2 μ L	3 μ L	4 μ L	5 μ L
0.5mL	10 μ L	20 μ L	30 μ L	40 μ L	50 μ L
1mL	20 μ L	40 μ L	60 μ L	80 μ L	100 μ L
1.25mL	25 μ L	50 μ L	75 μ L	100 μ L	125 μ L
2.5mL	50 μ L	100 μ L	150 μ L	200 μ L	250 μ L
5mL	100 μ L	200 μ L	300 μ L	400 μ L	500 μ L
12.5mL	250 μ L	500 μ L	750 μ L	1000 μ L	1250 μ L
25mL	500 μ L	1000 μ L	1500 μ L	2000 μ L	2500 μ L
50mL	1000 μ L	2000 μ L	3000 μ L	4000 μ L	5000 μ L

